

SteriPro 22B Vacuum Autoclave
Operator Manual - Model SP-22B

1. INTENDED USE

The SteriPro 22B is a class B vacuum steam steriliser for dental instruments. It has a 22 litre chamber and four validated cycles.

2. SAFETY

WARNING: The chamber door interlock must never be defeated. The chamber is a pressure vessel and opening it under pressure will cause serious scalding injury.

WARNING: Isolate the unit at the wall before removing any panel. Servicing behind the rear panel is for a qualified engineer only.

3. DAILY OPERATION

- 3.1 Fill the clean water reservoir to the MAX line with distilled water only. Tap water will scale the heater and void the warranty.
- 3.2 Empty the waste reservoir at the start of every session.
- 3.3 Load instruments in the trays supplied, hinged instruments open, with a 10 mm gap between pouches so steam can circulate.
- 3.4 Select the cycle: 134 C Universal for wrapped loads, 134 C Fast for unwrapped solid instruments, 121 C Delicate for handpieces.
- 3.5 Confirm the printout shows the plateau held for the full sterilisation hold time before releasing the load.

4. TROUBLESHOOTING

E01 - Water low. Refill the clean reservoir and restart the cycle.

E02 - Waste full. Drain the waste reservoir using the front tap.

E04 - Door not sealed. Reseat the load, check nothing fouls the gasket, wipe the gasket with a damp lint free cloth and retry.

E07 - Vacuum test failed. Run the leak test with an empty chamber after a cold start. If it fails twice, call the service engineer.

E11 - Temperature sensor fault. Stop using the unit and call the service engineer. Do not attempt to recalibrate the sensor.

5. SERVICE

Annual pressure vessel inspection and a service every 12 months or 1000 cycles, whichever is sooner. Record every service in the practice equipment log.